



INFORMATION WAVE SYSTEM

AVATAR

USER MANUAL

A personal information-wave system based on injection generators made from indium phosphide (InP) heterostructures. Concept, modes, applications and a recommended usage routine.

Contents

01 Device Concept

02 What's in the Box

03 Operating Modes

04 Working with Biomaterials

05 Using the Personalised Tablet

06 Water Structuring

07 Night Programme

08 Recommended Usage Routine

09 Safety Precautions

01

Device Concept

AVATAR is a personal information-wave system based on injection generators made from indium phosphide (InP) heterostructures.

The device is designed for:

- building the user's individual information profile;
- working with biomaterials;
- structuring water;
- local action on biologically active zones;
- supporting the body's natural self-regulation processes.

HOW IT WORKS

The technology is based on low-intensity electromagnetic processes and special semiconductor structures.



02

What's in the Box

01

AVATAR main unit

02

"Tablet" (×1)

03

Charging cable

04

User manual

INSIDE THE "TABLET"

The tablet contains five injection generators based on indium phosphide.



03

Operating Modes

MODE 1

Recording & clearing
information

1 min

DURATION

- preparing the tablet;
- recording information;
- structuring water.

MODE 2

Quick local action

5 min

DURATION

- during walks;
- when worn on the body;
- when working with biologically active zones.

MODE 3

Deep work

10 min

DURATION

- for prolonged wear;
- night programmes;
- working with individual information.



04

Working with Biomaterials

4.1 Recording information from saliva

Morning saliva, collected right after waking, is most commonly used.

Procedure

- 01 Place a small amount of saliva into an airtight bag.
- 02 Put the bag on the reading surface of the tablet.
- 03 Start the recording mode.
- 04 Wait for the cycle to finish.
- 05 Remove the bag.

RESULT

After recording, the tablet is considered personalised.

4.2 Recording information from urine

The first morning urine is used.

Procedure

- 01 Place a small amount of urine into a clean plastic cup.
- 02 Place the cup on the tablet.
- 03 Start the recording mode.
- 04 Wait for the cycle to finish.
- 05 Remove the sample.

RESULT

Once recording is complete, the tablet holds a personal information profile of the body's current state.

4.3 Working with sweat

For athletes and people engaged in intense physical activity.

You can:

- use a wipe with sweat;
- place it into a bag;
- record it the same way as saliva.



05

Using the Personalised Tablet

After recording, the tablet can be used in several ways.

5.1 Wearing on the body

The tablet can be placed:

- in a pocket;
- on a belt;
- in a chest pocket;
- in a dedicated case.

6–12 HOURS recommended time

5.2 During a walk

One of the most popular methods.

Recommended:

- record the morning information;
- take the tablet with you;
- take a 30–90 minute walk.

5.3 Targeting specific areas

The tablet can be applied to the area of the:

- kidneys;
- liver;
- solar plexus;
- spine;
- joints.

20–60 MINUTES per area

Water Structuring

OPTION 1**Basic structuring**

- 01 Use a cleared tablet.
- 02 Place a bottle of water on the tablet.
- 03 Start mode 1.
- 04 After one minute the water is ready.

This mode is used for everyday application.

OPTION 2**Personalised water**

- 01 First perform the basic structuring.
- 02 Then place the bottle on the personalised tablet.
- 03 Leave it for 6–8 hours.

In the morning the water is ready to drink.



07

Night Programme

In the evening:

- 01 Record from saliva.
- 02 Place the tablet under the sheet — over the area of the body where the organs you have complaints about are located and which need a wellness effect.
- 03 Leave it overnight.

IN THE MORNING

You can use the same tablet to structure water.



Recommended Usage Routine

TIME	ACTIONS
Morning	• record from saliva; • structure water; • drink water on an empty stomach.
Daytime	• wear the tablet; • a 30–60 minute walk.
Evening	• place locally on the area of interest; • prepare for the night programme.



Safety Precautions

IMPORTANT

- Avoid direct contact of liquids with the device.
- Use only airtight bags or containers.
- Do not immerse the device in water.
- Do not subject it to impacts or overheating.
- Store in a dry place.

